

# FUJI DRI-CHEM REFERENCE FLUID RE

## [Warnings and precautions]

1. Stop using when the volume of the reference fluid has decreased and meniscus level has gone below the lowest limit line indicated on the label or if three months have passed after first opening.
2. Seal up tightly after use to avoid concentration.
3. Warm up to room temperature before use.
4. In case of ingestion, eye or skin contact, wash with plenty of water immediately and consult medical personnel.

## [Intended use]

This product is used as reference solution for testing Na, K and Cl on FUJI DRI-CHEM system.

## [Limitations of use]

This product is specially designed for FUJI DRI-CHEM system. This product cannot be applied to other systems.

## [Composition]

Sodium chloride	100mmol/L
Sodium bicarbonate	30mmol/L
Potassium dihydrogenphosphate	4mmol/L
Sodium azide (preservative)	0.01%
Polyvinylpyrrolidone	
Glycerin	

## [Procedure]

Consult "INSTRUCTION MANUAL" for FUJI DRI-CHEM ANALYZER.

## [Disposal information]

As this product contains 0.01% sodium azide as a preservative, dilute the residual liquid with plenty of water when disposing.

## [Storage and handling]

Storage condition: Between 2-8°C(35.6-46.4° F)

## [Contents]

8 mL × 6 bottles



Consult "INSTRUCTION MANUAL" for FUJI DRI-CHEM ANALYZER or at <http://www.fujifilm.com/products/medical/>



FUJIFILM Europe GmbH  
Heesenstrasse 31, 40549 Düsseldorf, GERMANY



FUJIFILM Corporation  
26-30, Nishiazabu 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN

