



Date of issue: 1/Feb/2013

FUJI DRI-CHEM DILUENT DL(CRP)

1904000

[Warnings and precautions]

Warm up to room temperature before use.

[Intended use]

This product is used for diluting CRP test samples applied to FUJI DRI-CHEM system.

[Limitations of use]

This product is specially designed for FUJI DRI-CHEM system. This product cannot be applied to other systems.

[Composition]

MES buffer solution 0.42% (MES: 2-morpholinoethanesulfonic acid)

Sodium azide 0.02%

[Procedure]

Consult "INSTRUCTION MANUAL" for FUJI DRI-CHEM ANALYZER.

[Disposal information]

As this product contains 0.02% sodium azide as a preservative, dilute the residual liquid with plenty of water when disposing.

[Storage and handling]

Storage condition: Between 2-8°C

[Contents]

32mLx1

Lot number



Manufacturer

Authorized representative in the European Community

In vitro diagnostic medical devices



Consult instructions for use



Use by



Consult "INSTRUCTION MANUAL" for FUJI DRI-CHEM ANALYZER or at <http://www.fujifilm.com/products/medical/>

FUJIFILM Europe GmbH
Heesenstrasse 31, 40549 Düsseldorf, GERMANY



FUJIFILM Corporation
26-30, Nishiazabu 2-Chome, Minato-ku, Tokyo 106-8620, JAPAN